

Human Recombinant OX1 Orexin Receptor Stable Cell Line

Cat. No. M00224

Version 06092014

I	INTRODUCTION	1
II	BACKGROUND.....	1
III	REPRESENTATIVE DATA.....	2
IV	THAWING AND SUBCULTURING.....	2
V	REFERENCES	3
	Limited Use License Agreement.....	4

I. INTRODUCTION

Catalog Number: M00224

Cell Line Name: CHO-K1/OX1

Gene Synonyms: OX1, HRTR1

Expressed Gene: Genbank Accession Number NM_001525; no expressed tags

Host Cell: CHO-K1

Quantity: Two vials of frozen cells (3×10^6 per vial)

Stability: 16 passages

Application: Functional assay for OX1 receptor

Freeze Medium: 45% culture medium, 45% FBS, 10% DMSO

Complete Growth Medium: Ham's F12, 10% FBS

Culture Medium: Ham's F12, 10% FBS, 400 μ g/ml G418

Mycoplasma Status: Negative

Storage: Liquid nitrogen immediately upon delivery

II. BACKGROUND

The orexin/hypocretin peptides orexin A and orexin B (also known as hcrt-1 and hcrt-2) are 33- and 28-amino acid peptides, respectively. They are preferentially expressed in hypothalamus. The orexins have a range of physiological functions including feeding control, energy metabolism, neuroendocrine function modulation, and sleep-wake cycle regulation. The two orexin receptor subtypes OX1 and OX2 both mediate the action of orexin-A and orexin-B and couple efficiently through Gq/11 to activate phospholipase C and lead to elevation of intracellular calcium.

§: GenScript employs a PCR-based method to test the mycoplasma. The test covers 11 of the most common strains of mycoplasma, (covering approximately 95% of *M. fermentans*, *M. hyorhinis*, *M. arginini*, *M. orale*, *M. salivarium*, *M. hominis*, *M. pulmonis*, *M. arthritidis*, *M. neurolyticum*, *M. hyopneumoniae* and *M. capricolum*) and one species *Ureaplasma* (*U. urealyticum*), with sufficient sensitivity and specificity.

III. REPRESENTATIVE DATA

Concentration-dependent stimulation of intracellular calcium mobilization by Orexin A in CHO-K1/OX1 and CHO-K1 cells




Figure 1. Orexin A-induced concentration-dependent stimulation of intracellular calcium mobilization in CHO-K1/OX1 and CHO-K1 cells. The cells were loaded with Calcium-4 prior to stimulation with a OX1 receptor agonist, Orexin A. The intracellular calcium change was measured by FlexStation. The relative fluorescent units (RFU) were plotted against the log of the cumulative doses (10-fold dilution) of Orexin A (Mean \pm SD, n = 2). The EC_{50} of Orexin A on OX1 in CHO-K1 cells was 22 nM. The S/B of Orexin A on OX1 in CHO-K1 cells was 30.

Notes:

1. EC_{50} value is calculated with four parameter logistic equation:

$$Y = \text{Bottom} + (\text{Top} - \text{Bottom}) / (1 + 10^{-(\text{Log}EC_{50} - X) \cdot \text{HillSlope}})$$
 X is the logarithm of concentration. Y is the response
 Y is RFU and starts at Bottom and goes to Top with a sigmoid shape.
2. Signal to background Ratio (S/B) = Top/Bottom

IV. THAWING AND SUBCULTURING

Thawing Protocol

1. Remove the vial from liquid nitrogen tank and thaw cells quickly in a 37°C water-bath.
2. Just before the cells are completely thawed, decontaminate the outside of the vial with 70% ethanol and transfer the cells to a 15 ml centrifuge tube containing 9 ml of complete growth medium.
3. Pellet cells by centrifugation at 200 x g force for 5 min, and remove the medium.
4. Resuspend the cells in complete growth medium.
5. Transfer the cell suspension to a 10 cm dish with 10 ml of complete growth medium.

6. Grow the cells in incubator with 37°C, 5 %CO₂.
7. Add antibiotic in the following day.

Sub-culturing Protocol

1. Remove the culture medium from cells.
2. Wash cells with PBS (pH=7.4) to remove all traces of serum that contains trypsin inhibitor.
3. Add 2.0 ml of 0.05% (w/v) Trypsin- EDTA (GIBCO, Cat No. 25300) solution into 10 cm dish and observe the cells under an inverted microscope until cell layer is dispersed (usually within 3 to 5 minutes).

Note: To avoid cells clumping, do not agitate the cells by hitting or shaking the dish while waiting for the cells detach. If cells are difficult to detach, please place the dish in 37°C incubator for ~2 min.

4. Add 6.0 to 8.0 ml of complete growth medium into dish and aspirate cells by gently pipetting.
5. Centrifuge the cells at 200 x g force for 5min, and remove the medium.
6. Resuspend the cells in culture medium and add the cells suspension to new culture dish.
7. Grow the cells in incubator with 37°C, 5 %CO₂.

Subcultivation Ratio: 1:3 to 1:8 weekly.

Medium Renewal: Every 2 to 3 days

V. REFERENCES

1. Esmaeil Akbari, Nasser Naghdi and Fereshteh Motamedi (2006) Functional inactivation of orexin 1 receptors in CA1 region impairs acquisition, consolidation and retrieval in Morris water maze task. *Behavioural Brain Research* 173(1): 47-52
2. Kukkonen JP et al. (2002) Functions of the orexinergic/hypocretinergic system *Am. J. Physiol.* 283: C1567-C1591
3. Sakurai, T., Amemiya, A., Ishii, M., Matsuzaki, I., Chemelli, R.M., Tanaka, H., Williams, S.C., Richardson, J.A., Kozloski, G.P., Wilson, S., Arch, J.R.S., Buckingham, R.C., Haynes, A.C., Carr, S.A., Annan, R.S., McNulty, D.E., Liu, W.S., Terrett, J.A., Elshourbagy, N.A., Bergsma, D.J. and Yangisawa, M. (1998) Orexins and orexin receptors: a family of hypothalamic neuropeptides and G-protein coupled receptors that regulate feeding behaviour. *Cell* 92: 573-585

GenScript USA Inc,

860 Centennial Ave.

Piscataway, NJ 08854

Toll-Free: 1-877-436-7274

Tel: 1-732-885-9188, Fax: 1-732-210-0262

Email: product@genscript.com

Web: <http://www.genscript.com>

For Research Use Only.

Limited Use License Agreement

This is a legal agreement between you (Licensee) and GenScript USA Inc. governing use of GenScript's stable cell line products and protocols provided to licensee. By purchasing and using the stable cell line, the buyer agrees to comply with the following terms and conditions of this label license and recognizes and agrees to such restrictions:

- 1) The products are not transferable and will be used at the site where they were purchased. Transfer to another site owned by buyer will be permitted only upon written request by buyer followed by subsequent written approval by GenScript.
- 2) The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party.
- 3) The products sold by GenScript are for laboratory and animal research purposes only. The products are not to be used on humans, for consumption, or for any unlawful uses.

GenScript USA Inc. will not assert against the buyer a claim of infringement of patents owned or controlled by GenScript USA Inc. and claiming this product based upon the manufacture, use or sale of a clinical diagnostic, therapeutic and vaccine, or prophylactic product developed in research by the buyer in which this product or its components has been employed, provided that neither this product nor any of its components was used in the manufacture of such product. For information on the use of this product for other purposes, contact Marketing Department, GenScript USA Inc., 120 Centennial Avenue, Piscataway, New Jersey 08840, U.S.A. Phone: 1-732-885-9188. Fax: 1-732-210-0262. Email: marketing@genscript.com.