



Overexpress CRLF2 and IL-7R

overexpress Interleukin 7 receptor and CRFL2 simultaneously in kit225 cells



Relative level of IL-7R in Kit225 cells before and post lentivirus infection



Lane 1: positive control (pLVX-Puro-CRLF2-IRES-IL-7R plasmid)
Lane 2: Kit225 cells
Lane 3: lenti-Kit225 cells
Lane 4: positive control (pLVX-Puro-CRLF2-IRES-IL-7R plasmid)
Lane 5: Kit225 cells
Lane 6: Lenti-Kit 225 cells