

Sample Name : α -Factor Mating Pheromone, yeast
 Cat. No. : RP01002
 Lot. No. : PE8851603

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow: 1 ml/min
 Wavelength: 220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
30.00	Pumps	Pump A B.Conc	95	
30.01	Pumps	Pump A B.Conc	5	
40.00	Pumps	Pump A B.Conc	5	
40.01	Controller	Stop		

<<Column Performance>>

<检测器 A>

Column : AlltimaTM C18 4.6 x 250 mm

Chromatogram



1 Det. A Ch1 / 220nm

Peak Table

检测器 A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	18.155	336594	50702	1.579
2	18.696	278181	47613	1.305
3	18.856	20566395	1843093	96.454
4	19.175	76724	13095	0.360
5	19.494	64675	8837	0.303
Total		21322569	1963340	100.000