


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 07/23/18
Time Processed : 17:28:20
Injection Volume : 0.2
Sample Name : Urocortin III, mouse
Cat. No. : RP10101
Lot. No. : PE1571807
Theoretical MW : 4172.97
Observed MW : 4173.00

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH