


Mass Spectrum



Sample Information

User : Gary
Month-Day Processed : 10/24/18
Time Processed : 14:19:06
Injection Volume : 0.2
Sample Name : Parathyroid Hormone (PTH) (1-34), Human
Cat No. : RP01001
Lot No. : PE0491810
Theoretical MW : 4117.74
Observed MW : 4117.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Equipment : GK11010007
Interface Bias : +4.5kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2ml/min
B.conc : 50%H2O/50%MeOH