


Mass Spectrum



Sample Information

User :	Gary	Equipment	: GK11010007
Month-Day Processed :	09/26/19	Interface Bias	: +4.5kV
Time Processed :	05:40:31 PM	Drying Gas Flow	: 5 L/min
Injection Volume :	0.3	T.Flow	: 0.2ml/min
Sample Name :	Thyrotropin-Releasing Hormone (TRH), Free Acid	B.conc	: 50% H <sub>2</sub> O / 50% MeOH
Cat No. :	RP10114		
Lot No. :	E4689919I		
Theoretical MW :	363.37		
Observed MW :	363.1		