


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 03/27/18
 Time Processed : 12:29:25
 Injection Volume : 0.4
 Sample Name : β -Amyloid (1-42), rat
 Cat No. : RP10013
 Lot No. : PE3081803
 Theoretical MW : 4418.10
 Observed MW : 4417.60

Interface

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50%H₂O/50%MeOH