

Sample Name: Parathyroid Hormone (PTH) (1-34), Human

Cat. No. : RP01001

Lot. No. : E0334319F

Time Processed : 9:28:27

Month-Day-Year Processed : 06/19/2019

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow: 1 ml/min

Wavelength: 220 nm


Time	Unit	Command	Value	Comment
0.01	Pumps	Pump A B.Conc	5	
25.00	Pumps	Pump A B.Conc	65	
25.01	Pumps	Pump A B.Conc	95	
31.00	Pumps	Pump A B.Conc	95	
31.01	Pumps	Pump A B.Conc	5	
40.00	Pumps	Pump A B.Conc	5	
40.01	Controller	Stop		

<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: G09010878



1 Det.A Ch1 / 220nm

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	13.613	126631	12341	0.694
2	14.283	151346	16138	0.829
3	14.427	17512433	1949577	95.970
4	14.650	457413	33608	2.507
Total		18247824	2011663	100.000