


# Mass Spectrum



## Sample Information

Acquired by : Gary  
Month-Day Processed : 04/24/19  
Time Processed : 9:08:17  
Injection Volume : 0.4  
Sample Name : Urocortin II, human  
Cat No. : RP10098  
Lot No. : PE7031901  
Theoretical MW : 4450.28  
Observed MW : 4450.0

Interface :ESI  
Nebulizing Gas Flow :1.5L/min  
CDL Temp :250°C  
Block Temp :200°C

Equipment: SS-CM-0273  
Interface Bias: +4.5 kV  
Drying Gas Flow: 5 L/min  
T.Flow: 0.2 ml/min  
B.conc: 50%H2O/50%MeOH