


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 11/03/17
Time Processed : 16:43:57
Injection Volume : 0.2
Sample Name : Flagelin22
Cat. No. : RP19986
Lot. No. : PE4471710
Theoretical MW : 2272.50
Observed MW : 2272.50

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH