


Mass Spectrum



Sample Information

Acquired by : Gary
 Month-Day Processed : 12/25/17
 Time Processed : 14:38:09
 Injection Volume : 1
 Sample Name : β-Amyloid (42-1)
 Cat No. : RP10034
 Lot No. : PE6851710
 Theoretical MW : 4514.40
 Observed MW : 4513.80

Interface : ESI
 Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250°C
 Block Temp : 200°C

Interface Bias
 Drying Gas Flow
 T.Flow
 B.conc

: +4.5 kV
 : 5 L/min
 : 0.2 ml/min
 : 50%H₂O/50%MeOH