

Sample Name :Thyrotropin-Releasing Hormone (TRH)  
 Cat. No. : RP10112  
 Lot. No. :E4689519I  
 Time Processed :8:22:32  
 Month-Day-Year Processed :09/18/2019

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
 Total Flow:1 ml/min  
 Wavelength:220 nm

Time	Unit	Command	Value	Comment
0.01	Pumps	Pump B Conc.	5	
25.00	Pumps	Pump B Conc.	65	
25.01	Pumps	Pump B Conc.	95	
27.00	Pumps	Pump B Conc.	95	
27.01	Pumps	Pump B Conc.	5	
35.00	Pumps	Pump B Conc.	5	
35.01	Controller	Stop		


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment :SS-CM-0310

Chromatogram



1 Det.A Ch1 / 220nm

Peak Table

Detector A Ch1 220nm

Peak#	Ret. Time	Area	Height	Area %
1	5.489	39204014	1732879	99.725
2	7.708	39405	1354	0.100
3	8.519	44852	4927	0.114
4	8.689	9309	1683	0.024
5	9.329	8837	992	0.022
6	20.573	5534	617	0.014
Total		39311950	1742452	100.000