


# Mass Spectrum



## Sample Information

User : Gary  
Month-Day Processed : 01/31/20  
Time Processed : 06:27:25 PM  
Injection Volume : 0.3  
Sample Name : Urocortin II, human  
Cat No. : RP10098  
Lot No. : E5492920A  
Theoretical MW : 4450.28  
Observed MW : 4450.8

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : GK11010007  
Interface Bias : +4.5kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2ml/min  
B.conc : 50%H<sub>2</sub>O/50%MeOH