


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 12/25/17
Time Processed : 14:38:09
Injection Volume : 1
Sample Name : ?2-Amyloid (42-1)
Sample ID : U9173CJ160-1
Theoretical MW : 4514.40
Observed MW : 4513.80

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250°C
Block Temp : 200°C

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H2O/50% MeOH